Graphic Resources

Lab 1 (Fingerprints) Sc9 and 10

```
For each room:
```

6 Inkpads

12 dusting brushes

12 sellotapes

Paper towels

Microscopes slides

Corn starch (6 pots)


Talcum powder (6 pots)

Charcoal (6 pots)

12 laminated fingerprint pattern sheets

6 doorknobs with thumb print.

Hand lenses x 12













Right Thumb



SUSPECT A










Right Thumb





SUSPECT F












SUSPECT ****



























Lab 2 (Hair Analysis)

12 microscopes


Slides labelled 1-3 (6 sets)


Microscope slides


Sellotape x 6


Laminated identification sheets x 6


Animal Hair identification Sheet



















Animal Hair identification Sheet Electron Microscope







Image source: Safterstein Criminalists, 5th Ed. Prentice Hall, 1995



Sample 1



Sample 2



Sample 3



Remember the door code is 1657 577.

Lab 3 (Document Analysis)

| Suspect | Sample |
|---------|----------------------------------|
| Α | Remember the door code is 1657 |
| В | Remember the door code is 1657 |
| С | Remember the door code is 1657 |
| D | Remember the devor code is 1657 |
| E | Remember the door code is 1657. |
| F | Remember the dow code is 1657. |
| G | Remember The closer code is 1657 |
| Н | Remember the door code is 1657 |

Ava

Clair

Claire


Gordon

Gemma


lan

Jon M

Victoria




LF57




Forensic Crime Scene Science

LB62



LB72



Lab 4 (Ballistic Analysis)

Equipment required:

Examples of clothing from the suspects

Each group:

Sample tubes labelled A-H and each containing cotton buds.

(dirty the buds by placing in a plastic bag with a little charcoal and shaking)

Positive samples need to be impregnated with Potassium Nitrate: A, C, D, E, F and G

Negative samples: B and H

Beakers to place each sample set in.

Lead Nitrate dropper bottles (one per group)

Jigsaw puzzles of each casing

Laminated picture of casing found at the crime scene

